

A new development in control technology, grip chips enable users to easily activate quick grips, my grips or favorites for **i-limb ultra revolution**. The intuitive solution combined with the most dexterous hand available enables users to perform a variety of activities of daily living with ease and precise control.


### Features:

- New control technology offering more choices and enhanced ability to use **i-limb ultra revolution**<sup>1</sup>
- Placing the **i-limb ultra revolution** device in vicinity of grip chip activates grip through Bluetooth® connectivity
- my grips allows users to create 12 custom grips<sup>2</sup>
- Programming and changes made within the **my i-limb** or **biosim** apps
- Easy setup and minimal training required



### How grip chips work:

Using the **my i-limb** app, patients can easily program grip chips with any grip or gesture and place in a desired location. Simply holding an **i-limb ultra revolution** device within close range to a grip chip will activate the programmed feature.



Product Number	Description
PL238089	grip chip Starter Kit (1-10)
PL238090	grip chip Replenishment Kit (11-20)



<sup>1</sup>Additional software purchase and activation is required to use grip chips.

<sup>2</sup>For more information on my grips, reference the my i-limb app documents online at: <http://www.touchbionics.com/downloads/document-library>.